

NIV face mask (Non-vented Full Face Mask)

Good Vision

Has no forehead support and does not interfere with vision

Easy to Use

Easily adjusted and comfortable headgear

Improved Comfort

The soft High- Quality silicone mask offers minimal facial contact for a full-face mask and does not irritate the skin


Safe Material

Contains no PVC, DEHP, or LATEX




SPECIFICATION

PRODUCT PARTS



- Materials: silicone mask; polycarbonate frame
- Monitoring port: 4.3mmOD, meets ISO 594-1

HEADGEAR ASSEMBLY DIAGRAM



- L-type connector: 22mmID, meets ISO 5356-1
- Headgear: elastic fiber

ORDER INFORMATION

REF	Description	Size	Package
1701	NIV Face Mask	L	1set/box, 20sets/carton
1702	NIV Face Mask	M	1set/box, 20sets/carton
1703	NIV Face Mask	S	1set/box, 20sets/carton



DA CHUNG MEDICAL CO., LTD.

No.600, 602, Sec.1, Dong Shan Rd., Dong Shan Township, Yilan County 26948, Taiwan

TEL : +886 3 959 4777

FAX : +886 3 959 4707

E-mail : info@dc-med.com.tw

