

Product Name: Oxygen Tubing

Function :

The airway connector is a device used to connect a source of breathing gas to an endotracheal tube, a tracheotomy tube, or a mask.

Instruction for Use :

1. Connect the oxygen tube connector to the intended equipment and verify that the oxygen tube connector is indeed engaged with the equipment.
2. Use according to the intended use of the equipment (and follow the instructions of the individual equipment instructions).

Warning :

1. The product is used by a single patient.
2. Please make sure all parts are connected before use.
3. This instruction is only a general indication, and it still needs to be used by professional medical personnel.
4. Please confirm the appearance of the package and the product before use. If it is damaged or damaged, please do not use it.
5. The maximum recommended life of this device is 7 days.
6. Dispose of in accordance with the local medical waste disposal regulations.

Specification : Blue/Green 、 2M

Valid period: 3 years



Obelis s.a
Bd. Général Wahis 53 1030 Brussels, BELGIUM



Manufacturer : DA CHUNG MEDICAL CO., LTD.
Address : No. 44, Te-Xing 6th Rd, Dong-Shan
Township, Yi-Lan 26950 Taiwan.
Website : <http://www.dc-med.com.tw>

REV 01