

CPAP comfort +



YOUR BEST PARTNER



D.C. Medical
Solution for your sleep dyspnea

Make sleep joyful again
with
D.C. Medical CPAP comfort +



DA CHUNG MEDICAL CO., LTD.

No.600, Sec.1, Dong Shan Rd., Dong Shan Township, Yilan County 26948, Taiwan
TEL: +886 3 959 4777 FAX: +886 3 959 4707 E-mail: info@dc-med.com.tw



More than 80% of obstructive sleep apnea patients are not yet diagnosed¹...

Do you have²...

- Daytime sleepiness
- Snoring
- Pause in breathing during sleep
- Difficulties with memory and concentration
- Unusual moodiness or irritability
- Frequently waking up to urinate at night
- Morning headaches
- Dry mouth

Treating sleep apnea for better health now⁴

Research has found association between sleep apnea and...

- Type 2 diabetes
- Stroke
- Heart attack
- Shortened lifespan

“CPAP is the first-line treatment, and it’s very effective.”

By Jonathan Jun, M.D. , a pulmonary and sleep medicine specialist at the Johns Hopkins Sleep Disorders Center.

Our CPAP comfort + can be the solution³

1. Sleep Sci. 2021 Jul-Sep; 14(3): 257–265.

2. American Lung Association Scientific and Medical Editorial Review Panel.

3. Product IFU, IFU-1601 REV 03

4. The Johns Hopkins University, The Johns Hopkins Hospital, and Johns Hopkins Health System. Article: The Dangers of Uncontrolled Sleep Apnea (<https://www.hopkinsmedicine.org/health/wellness-and-prevention/the-dangers-of-uncontrolled-sleep-apnea>) last assessed 17-Arp-2023

Never easy like this, total facial pressure is 12-38% less than other generics⁵

Sturdy, high elasticity and breathable mesh with easy on-off velcro design.

Study proven bio-compatible fabric, zero irritation, zero allergy.

Fully silicone made, provide you top notch comfort experience.

PC hard shell, impact durable, long lifespan

Quick buckle design, zero interference to diet and night activities.

Solid material, non flake, zero inhalation risk.

Unique auto-fit design, facial size and shape wouldn't get in the way


No metallic component, MRI ready, safe for pacemaker and other electronic devices

Anti-suffocation valve, safety assured


360 degrees rotatable connector, provide nature sleep experience



PRODUCT PARTS



HEADGEAR ASSEMBLY DIAGRAM



- Materials: silicone mask; polycarbonate frame
- L-type connector: 22mmOD, meets ISO 5356-1
- Monitoring port: 4.3mmOD, meets ISO 594-1
- Headgear: elastic fiber

REF	Description	Size	Package
1631	CPAP Face Mask (double layers)	L	1set,box, 20sets/carton
1632	CPAP Face Mask (double layers)	M	1set,box, 20sets/carton
1633	CPAP Face Mask (double layers)	S	1set,box, 20sets/carton



Find your size to resume quality sleep

Instructions for use:

Step 1: Align the red dot with your bottom lip, as shown in the picture.

Step 2: Place edge of the card between your eyes to decide the suitable mask size.

med

YOUR BEST PARTNER

REF	Product		
1601 / 1701	CPAP / NIV FACE MASK - LARGE		LARGE
1602 / 1702	CPAP / NIV FACE MASK - MEDIUM		MEDIUM
1603 / 1703	CPAP / NIV FACE MASK - SMALL		SMALL

CTN-01-109-01 VER.01